

SALVO+GEL

INSTANT HAND SANITIZER

Germs remain active on hard surfaces for hours, and sometimes even days. These germs are passed on from one person to another when we handle everyday things like telephones, computer keyboards and mice, door handles, hand rails, money, etc. Experts acknowledge that approximately 80% of all common infections are spread by the hands this way.

Opus SALVO-GEL is an effective weapon in the ever more serious battle against pathogenic bacteria.

Opus SALVO-GEL has been developed to meet the requirements of the food industry for a completely effective, yet quick-acting and nontainting, food-safe hand steriliser. Combining a high level of alcohol gives proven bactericidal effectiveness.

Added glycerine leaves hands soft, supple, non-sticky and moisturised - does not dry out the skin and leaves no after-odour.

These advanced properties make Opus SALVO-GEL equally suitable for many medical and care service applications - in fact, any application where there is a risk of cross-infection from hand contact.

Independent laboratory tests prove 99.9999% kill of gram positive and gram negative bacteria.

Complies with EU Approval EN1276 for food use and is effective against *Pseudomonas Aeruginosa*, *Escherichia Coli* ATCC 10536, *Staphylococcus Aureus* ATCC 6538, and *Enterococcus hirae* after a contact time of 5 minutes.

Complies with EU Approval EN13727 for effective use against MRSA, *Salmonella Typhimurium*, *Listeria Monocytogenes* in medical areas.

Apply a small quantity to clean, dry hands and rub well in, paying particular attention to fingernails and creases in the skin. Within 15 to 20 seconds, SALVO-GEL will have evaporated, leaving the skin soft and dry. It may also be used in the same way to sterilise surgical gloves.

COSHH Data: An EC safety data sheet (MSDS) providing additional information is available on request.

Packaging Sizes:

- 24 x 80 ml pump or flip-top
- 12 x 250 ml pump or flip-top
- 10 x 500 ml pump bottles
- 6 x 1 litre pump bottles
- 6 x 1000 ml pouches for use in 1 litre dedicated dispenser
- 6 x 2000ml pouches for use in 2 litre pouch dispenser
- 4 x 5 litre bulk-fill jerrycans

SALVO-GEL IS USED IN:

- Clinics
- Dental surgeries
- Doctors' surgeries
- Food preparation areas
- Hospitals
- Kindergartens
- Nursing homes
- Schools
- Colleges
- Butchers' shops
- Veterinary surgeries
- Local Councils
- Supermarkets
- Bakeries



Excalibur Hygiene Limited - Opus Handcare Division
 Unit 42 Enterprise Trading Estate | Pedmore Road
 Brierley Hill | West Midlands | DY5 1TX | United Kingdom
 Tel: +44 (0)1384 671505 | Fax: +44 (0)1384 671256
 Email: help@excalibur-hygiene.co.uk
 Website: www.hand-cleaners.co.uk

